



## Application Instructions – Rebound Walking Boot

This product has been designed and tested based on single patient usage. It is a single-use device and should never be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

### Cautions/Warnings

READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE USING THE DEVICE.

- Federal law requires this device to be sold by a physician or by the order of a physician.
- Ensure not to over inflate air chamber.
- Caution should be used if prescribed for patients with diabetes, vascular deficiency, and neuropathy:
  - o The skin should be checked regularly
  - o Over inflation of air chamber may cause skin compromise in patients with diabetic neuropathy. Deflate air chamber to reduce compression with patients reporting discomfort or showing skin discoloration.
- Ensure patient performs circulation checks. If a loss of circulation is felt, or if the walker feels too tight, deflate air chamber to reduce compression to a comfortable level. If discomfort continues, the patient should contact their medical professional immediately.
- This device should not be applied on patients who are not able to communicate discomfort.
- In high altitude conditions, the air chamber will expand past their optimal level. Deflate air chamber to reduce compression to a comfortable level.
- User should use caution when walking on slippery or wet surfaces to avoid injury.
- User maximum weight limit 300 lbs.

### Preparation

1. Open liner by first unfastening straps and removing front shell.

Note: Do not remove air bulb from liner.

### Application

2. While seated, place leg in brace.
3. Wrap ankle pad snugly around ankle.  
Note: Ankle pad can be cut to size.
4. Wrap liner securely around foot then leg.
5. Place front shell over ankle making sure that the air bulb is aligned with the opening.
6. Fasten ankle, foot, and calf strap in this order. Tighten straps until snug and comfortable.

### Inflation/Deflation


7. Inflate air chambers by pressing "+" on the air bulb, and deflate by pressing "-". Inflate or deflate until snug. Caution: Do not overinflate!
8. Before removing brace, deflate air chambers by pressing "-".

### Washing Instructions

- Remove liner from walker before washing
- The soft good liner should be washed by hand using a mild detergent and cold water. Do not use bleach or strong detergent. Do not wash in washing machine or place in dryer.
- Air dry flat.
- Walker plastic can be washed with a damp cloth and mild soap.



  
reddot design award  
winner 2011



*Moisture-wicking softgood  
helps keep skin dry  
and comfortable*

*Linked medial and lateral air bladders  
with single inflation point  
provide balanced compression and are  
shaped to cradle the Achilles, holding the  
heel in place*

